



**Be Right™**



## **Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Blue Cap/ Stem, 6/pk**

**Product #:**

**2928202**

**ZAR Price:**

**Contact Hach**

**Ships within 1 week**

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

---

### **Specifications**

Capacity:	500 mL
Color:	Clear
Material:	LDPE
Quantity:	6 pk
Stem:	Blue Cap/Stem
Unit:	6 pk