



Be Right™



Gas washing bottle, 1200 mL

Product #:	2662200
ZAR Price:	Contact Hach
Call for ship date	

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis—Hach Method 10065.

Specifications

Body Material:	Glass
Capacity:	1200 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Quantity:	1 each