



Test Tube, rack, polyethylene, 20 mm (O.D.), 20 holes

| | |
|-------------------|----------------|
| Product #: | 2497912 |
| ZAR Price: | Contact Hach |
| Available | |

Twenty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may Vary.

Specifications

| | |
|------------------------|--------|
| Diameter: | 20 mm |
| Holes: | 20 |
| Quantity: | 1 each |
| Test Tube Arrangement: | 4 x 5 |