



Test tube/ MPN rack, polyethylene, 20 mm (O.D.), 40 holes

Product #:	2497902
ZAR Price:	Contact Hach
Ships within 2 weeks	

Forty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifications

Diameter:	20 mm
Holes:	40
Quantity:	1 each
Test Tube Arrangement:	4 x 10