



Test tube, rack, polyethylene, 16 mm (O.D.), 60 holes

Product #:	2497901
ZAR Price:	Contact Hach
Ships within 2 weeks	

Sixty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may vary.

Specifications

Diameter:	16 mm
Holes:	60
Quantity:	1 each
Test Tube Arrangement:	5 x 12