



## Rack, 21 x 1 inch round cuvettes

<b>Product #:</b>	<b>2497904</b>
ZAR Price:	Contact Hach
Available	

Twenty-one place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may vary.

---

### Specifications

Diameter:	30 mm
Holes:	21
Quantity:	1 each
Test Tube Arrangement:	3 x 7