



Rack, Test Tube, Polyethylene, 13 mm (O.D.), 90 Holes

Product #:	2497900
ZAR Price:	Contact Hach
Ships within 2 weeks	

Ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifications

Diameter:	13 mm
Holes:	90
Quantity:	1 each
Test Tube Arrangement:	6 x 15