



Pipet tips, TenSette pipet 1970010, 1.0-10.0 mL, non-sterile, pk/250

Product #:	2199725
ZAR Price:	Contact Hach
Ships within 1 week	

Perform accurate, economical pipetting with Hach's TenSette pipets. Single-channel, air-displacement pipets come in two models: 0.1-1.0 mL volumes (0.1 mL increments) or 1.0-10.0 mL (1.0 mL increments). Volume selection with a simple twist of the turret.

Specifications

Quantity:	250 /pk
-----------	---------