



Pipet, serological, sterile, polystyrene, 5 mL, blue, 200/pk

Product #:	2092667
ZAR Price:	Contact Hach
Call for ship date	

Sterile, disposable plastic pipet. Individually wrapped. Color coded. Calibrated “To Deliver”.

Graduations are calibrated for accurate dispensing to within ±2%

Pipets meet BST/TSE requirements; latex-free

DNase/RNase-free, non-pyrogenic, non-cytotoxic, non-haemolytic

Specifications

Capacity:	5 mL
Graduations:	0.10 mL
Negative Graduations:	2.0 mL
Quantity:	200 /pk