



Pipette, variable, volume 1.0-5.0 mL

Product #:	BBP065
ZAR Price:	Contact Hach
Ships within 1 week	

Hach's variable volume single-channel, air-displacement pipettes feature robust design, comfortable ergonomics, slim-line tip ejection, and single-handed operation. Outstanding accuracy with adjustments in 50 μ L increments (0.05 mL). Ideal for use with Hach LCK reagents. Supplied with calibration report and maintenance tool.

Specifications

Quantity:	1 each
What's included?:	Pipette 1.0-5.0 mL, maintenance tool, calibration certificate, basic manual.
	Tips are not included and have to be purchased separately.

What's included?

Pipette 1.0-5.0 mL, maintenance tool, calibration certificate, basic manual. Tips are not included and have to be purchased separately.