



Set of 2 pipettes, variable volume, incl. tips

Product #:	LZP320
ZAR Price:	Contact Hach
Ships within 2 weeks	

Value pack, consisting of two variable pipettes: BBP078 (0.1 - 1.0 mL) and BBP065 (1.0 - 5.0 mL), incl. tips

Specifications

Quantity:	1 each
Sample size:	0.1 - 1.0 mL / 1.0 - 5.0 mL
What's included?:	BBP078 Variable Pipette 0.1 - 1.0 mL
	BBP079 Pipette Tips for BBP078, 100/pk
	BBP065 Variable Pipette 1.0 - 5.0 mL
	BBP068 Pipette Tips for BBP065, 75/pk

What's included?

BBP078 Variable Pipette 0.1 - 1.0 mLBBP079 Pipette Tips for BBP078, 100/pkBBP065 Variable Pipette 1.0 - 5.0 mLBBP068 Pipette Tips for BBP065, 75/pk