



## Pipet, volumetric, class A, certified, 10 mL, 6/pk

|            |              |
|------------|--------------|
| Product #: | 2636168      |
| ZAR Price: | Contact Hach |
| Available  |              |

Glass Class A pipets serialized and calibrated “To Deliver”. Each comes with Certificate of Identification and Capacity. Marked with an “A” to designate compliance to ASTM E969.

---

### Specifications

|           |       |
|-----------|-------|
| Capacity: | 10 mL |
| Quantity: | 6 /pk |