



Pipet, Volumetric, Class A, Certified, 25 mL, 6/pk

| | |
|----------------------|----------------|
| Product #: | 2636170 |
| ZAR Price: | Contact Hach |
| Ships within 2 weeks | |

Glass Class A pipets serialized and calibrated “To Deliver”. Each comes with Certificate of Identification and Capacity. Marked with an “A” to designate compliance to ASTM E969.

Specifications

| | |
|-----------|-------|
| Capacity: | 25 mL |
| Quantity: | 6 /pk |