



## Bottle Carrier, Linear

<b>Product #:</b>	<b>2094300</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Linear, low-density polyethylene. Securely holds up to six standard 0.5 L glass bottles, pint-size reagent bottles, and most other 500 mL glass or plastic bottles—making breakage almost impossible. Lightweight, chemical-resistant, and easy to clean. Dimensions: 34.9 x 24.8 x 27.6 cm (14 x 10 x 11"). Maximum bottle size it will hold: 73 mm in diameter, 178 mm in height.

---

### Specifications

Depth:	27.6 cm
Length:	349 mm
Quantity:	1 each
Width 2:	10 Inches