



## Stainless Steel Replacement K-Type Probe

<b>Product #:</b>	<b>2630601</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Light weight and inexpensive, this K & J type thermometer provides extremely accurate and reliable temperature measurement. Features relative minimum and maximum temperature, hold, auto power off, and touch tone buttons. Includes one K-type probe. Probe Range (K): -328 to 2498°F Dimensions: 5 1/4" x 2 1/4" x 1 1/4" Weight: 6 oz.

---

### Specifications

Quantity:	1 each
Temperature Range:	-328 - 2498 °F