



Thermometer, Enviro-Safe, Non-Mercury, 0 to 230°F

Product #: 2635703

OBSOLETE ITEM
This item is no longer available.

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Includes a multi-language (English, French, German, Italian, Portuguese, and Spanish) Statement of Accuracy and Traceability

Specifications

Length:	305 mm
Quantity:	1 each
Temperature Range:	0 °F - 230 °F
Thermometer Divisions:	2 °F
Thermometer Fluid Type:	Non-Mercury