



Thermometer, NIST-Certified

Product #: 2551500

OBSOLETE ITEM
This item is no longer available.

Certified to NIST standards at specific temperature points, are supplied with a Calibration Report stating corrections to one-tenth of the smallest division and feature individual serial numbers, accuracy of $\pm 2\text{ }^{\circ}\text{C}$; $\pm 2.5\text{ }^{\circ}\text{C}$ above $110\text{ }^{\circ}\text{C}$; $\pm 3\text{ }^{\circ}\text{C}$ above $200\text{ }^{\circ}\text{C}$. Blue organic liquid against white back glass. Tested and certified from supplier's accredited/registered ISO/IEC 17025:2005, A2LA accredited laboratory against equipment whose calibration is traceable to NIST. Individually packaged NIST Certified.

Certified at 0, 10, 20, 30, 37, 40, and $50\text{ }^{\circ}\text{C}$.

Specifications

Length:	460 mm
Quantity:	1 each
Temperature Range:	$\pm 2\text{ }^{\circ}\text{C}$, $\pm 2.5\text{ }^{\circ}\text{C}$ above $110\text{ }^{\circ}\text{C}$, $\pm 3\text{ }^{\circ}\text{C}$ above $200\text{ }^{\circ}\text{C}$
Thermometer Divisions:	$0.1\text{ }^{\circ}\text{C}$
Thermometer Fluid Type:	Organic Filled