



## Thermometer, Enviro-Safe, Non-Mercury, -10 - 225 °C

<b>Product #:</b>	<b>2635700</b>
ZAR Price:	Contact Hach
Available	

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO 9000 or other quality systems.

---

### Specifications

Length:	405 mm
Quantity:	1 each
Temperature Range:	-10 °C - 225 °C
Thermometer Divisions:	1 °C
Thermometer Fluid Type:	Non-Mercury