



## Weight Set, Alloy 8 Metric Precision Calibration, ASTM Class 1, 100g - 1mg

Product #:

2617600

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

ASTM Class 1 Troemner Alloy 8 Highly polished stainless steel weight set. 100g-1g, 500 mg down to 50mg weights are 304 stainless steel, weights below 50mg are flat aluminum weights. Polypropylene case and forceps included. NVLAP Accredited Calibration Certificate included.

---

### Specifications

Quantity: 1 each