



Weight Set, Class M-2

Product #:	2101100
ZAR Price:	Contact Hach
Ships within 3-5 weeks	

Standard weight set. 16-piece, with stainless steel forceps. Eight fractional weights: 10, 20, 50, 100, 200, 200, 200, and 500 mg. (Two fractional weights, 200 mg and 500

mg, are stainless steel, the rest are aluminum) Eight brass gram weights: 1, 2, 2, 5, 10, 10, 20, and 50 grams. Polypropylene case.

Specifications

Quantity:	1 each
What's included?:	Includes: Brass weights and polypropylene case.

What's included?

Includes: Brass weights and polypropylene case.