



Field Balance, 2000 x 1.0g Capacity

Product #: 2946802

OBSOLETE ITEM
This item is no longer available.

Simple and easy to operate, this rugged field balance features a stability indicator, overload indicator, simple external calibration, large LCD display, auto power off, 2-year warranty, and four units of measure: g, oz, ct and tl. Includes convenient carrying case, stainless steel pan, AC power adaptor and 4 AA batteries. Overall dimensions: 5.1" x 7.6" x 2.0". Weight: 1.1 lb. Pan Size: 5.1" in diameter.

Compact

Small size means its easy to use anywhere. The scale can also be stored on its side or rear, saving space!

Simple and easy to operate

Three button operation and large LCD display make using the HLi series school, home or field friendly.

Multi weighing units

Meet any weighing application with weighing units of g, oz, lb-oz and tael for HL-2000i.

Accurate reading

Precision weighing of 2000g x 1g for HL2000i

Large weighing pan

Easily weigh large objects on the 130mm diameter weighing pan.

Specifications

Accuracy 2:	1 g
Capacity:	2000 g
Communication Capabilities:	None
Electrical Specifications:	AC power adaptor and 4 AA batteries
Length:	193 mm
Package size UOM:	each
Power:	4 AA Batteries
Quantity:	1 each
Resolution:	1.0 g
Width US:	5,1