



## Balance, Ohaus Model YA501

**Product #:** 2796900  
**ZAR Price:** Contact Hach  
**Ships within 1 week**

Designed for portability, this scale is ideal for remote or "on the go" applications. Weighs in g, oz, ct, gn, dwt. Unique hinged cover protects the scale with auto shut-off feature. Balance measures to 0.1 g. 500 g capacity. Weigh Pan measures 2.8 in x 2 in (70 mm x 50 mm)

### Protected pocket scale

Hinged cover is designed to cover the scale when not in use.

### Rugged design

Rugged plastic housing and stainless steel weighing pan provide a long-lasting portable scale.

### Units of measure

Weight units selectable from grams (g), ounces (oz) carats (ct) pennyweight (dwt) and grain (gn).

### Energy saving

Auto-off feature ensures longer battery life.

---

## Specifications

Accuracy 2:	0.1 g
Capacity:	500 g
Communication Capabilities:	None
Dimensions (H x W x D):	0.8 in x 2.9 in x 4.3 in (20 mm x 75 mm x 110 mm)
Power:	2 AAA Batteries
Quantity:	1 each
What's included?:	Includes portable scale and 2AAA Batteries.

---

## What's included?

Includes portable scale and 2AAA Batteries.