



## Field balance, 200 x 0.1 g capacity

**Product #:** 2946801

**OBSOLETE ITEM**  
This item is no longer available.

Simple and easy to operate, this rugged field balance features a stability indicator, overload indicator, simple external calibration, large LCD display, auto power off, 2-year warranty, and four units of measure: g, oz, ct and tl. Includes convenient carrying case, stainless steel pan, AC power adaptor and 4 AA batteries. Overall dimensions: 5.1" x 7.6" x 2.0". Weight: 1.1 lb. Pan Size: 5.1" in diameter.

### Compact

Small size means its easy to use anywhere. The scale can also be stored on its side or rear, saving space!

### Simple and easy to operate

Three button operation and large LCD display make using the HLi series school, home or field friendly.

### Multi weighing units

Meet any weighing application with weighing units of g, oz, ct and tael for HL-200i.

### Accurate reading

Precision weighing of 200g x 0.1g for HL-200i

### Large weighing pan

Easily weigh large objects on the 130 mm diameter weighing pan.

---

## Specifications

Accuracy 2:	0.1 g
Capacity:	200 g
Communication Capabilities:	None
Electrical Specifications:	AC power adaptor and 4 AA batteries
Length:	193 mm
Power:	4 AA Batteries
Quantity:	1 each
Resolution:	0.1 g
Width US:	5,1