



Lab Stirrer, 6-Paddle Jar Tester with 1 L Beakers

Product #: 2631710

OBSOLETE ITEM

This item is no longer available.

Phipps & Bird 900 Model Jar Tester allows you to replicate actual plant conditions and operations. Jar Tester runs in both standard and programmable modes. Features digital readout, stainless steel paddles placed 6 inches apart, regulated variable speeds of all paddles simultaneously, from 1 - 300 rpm, powder-coated steel uniframe chassis, built-in L.E.D. illuminator, anti-glare curtain, and dust cover. The unit holds either round 1000 mL mixing beakers or square 2 L Wagner Floc Jars.

Includes 1 L round glass beakers.

High accuracy and repeatability

PB Jar Testers utilize four operator-programmable memory banks

Quality craftsmanship

Built on a powder-coated steel uni-frame chassis

Versatile programming

Each memory bank can be run independently in any one of the memories, continuously, at a selected speed or automatic transitionin from a fast mix/short run to a slow mix/long run, or any combination or speeds and times desired to replicate actual plant conditions and operation

Long Lasting LED Illuminator

Provides reliable, diffused sample illumination

CE certified

Operates on 24VDC with a power supply that works with 100 - 240VAC electrical sources.

Specifications

Number of Paddles:	6
Options:	1L Round Beakers
Quantity:	1 each
Weight:	22.2 kg
What's included?:	Includes: 6 1L round glass beakers, dust cover, power supply and power cords.

What's included?

Includes: 6 1L round glass beakers, dust cover, power supply and power cords.