



## Lab Stirrer, 6-Paddle Jar Tester with 2 L Beakers

Product #: 2631720

**OBSOLETE ITEM**  
This item is no longer available.

Phipps & Bird 910 Model Jar Tester runs in both standard and programmable modes. Features digital readout, stainless steel paddles placed 6 inches apart, regulated variable speeds of all paddles simultaneously, from 1-300 rpm, powder-coated steel uniframe chassis, built-in L.E.D. illuminator, anti-glare curtain, and dust cover. The unit holds either round 1000 mL mixing beakers or square 2 L Wagner Floc Jars.

Includes 2L B-Ker<sup>2</sup> square acrylic jars.

### High accuracy and repeatability

PB Jar Testers utilize four operator-programmable memory banks

### Quality craftsmanship

Built on a powder-coated steel uni-frame chassis

### Versatile programming

Each memory bank can be run independently in any one of the memories, continuously, at a selected speed or automatic transitionin from a fast mix/short run to a slow mix/long run, or any combination or speeds and times desired to replicate actual plant conditions and operation

### Long Lasting LED Illuminator

Provides reliable, diffused sample illumination

### CE certified

Operates on 24VDC with a power supply that works with 100 to 240VAC electrical sources.

---

## Specifications

Number of Paddles:	6
Options:	2L Square Beakers
Quantity:	1 each
Weight:	22.2 kg
What's included?:	Includes: 6 2L square acrylic jars, dust cover, power supply and power cords.

---

## What's included?

Includes: 6 2L square acrylic jars, dust cover, power supply and power cords.