



Lab Stirrer, 4-Paddle Jar Tester with 2 L Beakers

Product #: 2631740
ZAR Price: Contact Hach
Ships within 3-5 weeks

Phipps & Bird 960 Model Jar Tester runs in both standard and programmable modes. Features digital readout, stainless steel paddles placed 6 inches apart, regulated variable speeds of all paddles simultaneously, from 1 - 300 rpm, powder-coated steel uniframe chassis, built-in L.E.D. illuminator, anti-glare curtain, and dust cover. The unit holds either round 1000 mL mixing beakers or square 2 L Wagner Floc Jars.

Includes 2 L B-Ker² square acrylic jars.

High accuracy and repeatability

PB Jar Testers utilize four operator-programmable memory banks

Quality craftsmanship

Built on a powder-coated steel uni-frame chassis

Versatile programming

Each memory bank can be run independently in any one of the memories, continuously, at a selected speed or automatic transition from a fast mix/short run to a slow mix/long run, or any combination of speeds and times desired to replicate actual plant conditions and operation

Long Lasting LED Illuminator

Provides reliable, diffused sample illumination

CE certified

Operates on 24VDC with a power supply that works with 100 to 240VAC electrical sources.

Specifications

Number of Paddles:	4
Options:	2L Square Beakers
Quantity:	1 each
Weight:	16.7 kg
What's included?:	Includes: 4 2L square acrylic jars, dust cover, power supply and power cords.

What's included?

Includes: 4 2L square acrylic jars, dust cover, power supply and power cords.