



## Oven, Laboratory, 120 Vac, 60 Hz

**Product #:** 1428900

**OBSOLETE ITEM**  
This item is no longer available.

For drying, baking, annealing, conditioning, sterilizing, evaporating, and dehydrating. Factory tested. Sturdy welded steel construction. Ample insulation. Two nickelplated shelves adjustable at 12.7 mm (0.5") increments. Thermostatically controlled with damper-controlled, gravity convection air circulation. Red pilot light indicates heater action. Gray-green, hammerloid, baked enamel finish. Three-wire cord and plug standard. Inside dimensions: 30w x 25d x 25h cm (11.8 x 9.8 x 9.8"). Outer dimensions: (19h x 14w x 12d"). Temperature range is ambient to 200°C. Control is  $\pm 15^\circ\text{C}$ .

### Compact

Small enough to fit on the lab bench.

### Economical

Perform general lab drying and baking at a low cost.

---

## Specifications

Accuracy 2:	$\pm 15^\circ\text{C}$
Dimensions (H x W x D):	483 mm x 356 mm x 305 mm
Electrical Specifications:	120 Vac, 50/60 Hz
Length:	483 mm
Quantity:	1 each
Weight:	Shipping: 18.5 kg
What's included?:	Included: Oven, three wired cord and plug, and two adjustable shelves.

---

## What's included?

Included: Oven, three wired cord and plug, and two adjustable shelves.