



Centrifuge, Laboratory, 8 x 15 mL 115 Vac, 60 Hz

Product #: 2676500

OBSOLETE ITEM
This item is no longer available.

Centrifuge supports small volume applications. The centrifuge is supplied with an angle rotor that can hold up to 8 x 15 mL, or 10 mL / 7 mL / 5 mL vacuum tubes with adapters (included).

Recommended for the centrifugation of samples for water, food, and soil analyses. This centrifuge is also used in universities and other research settings.

Versatile

Variable Speed: 300 - 4500 rpm \pm 20 rpm; increment: 100 rpm

Specifications

Electrical Specifications:	115 Vac, 60 Hz, 2.1 A
Power:	115 Vac, 60 Hz
Quantity:	1 each
Special Features:	Centrifuge, Laboratory
What's included?:	LCD Digital Clinical Centrifuge, with Rotor & Rotor Adapters, Adapter Plugs, power cord with US Plug, (110-220V, 50/60Hz).

Also included: Wrench for taking off the rotor, Philips screwdriver for manually opening the lid, operating manual, and QC certificate.

What's included?

LCD Digital Clinical Centrifuge, with Rotor & Rotor Adapters, Adapter Plugs, power cord with US Plug, (110-220V, 50/60Hz).Also included: Wrench for taking off the rotor, Philips screwdriver for manually opening the lid, operating manual, and QC certificate.