



## Stirrer, Mixer, Hand-Held

<b>Product #:</b>	<b>2496400</b>
ZAR Price:	Contact Hach
Available	

Hand-held mini-mixer used with test tubes, micro centrifuge tubes and containers too narrow for conventional mixers. Stirring rod flares out for use in larger spaces.

Clear vinyl carrying case, 2 AA batteries included.

Stirring rod dimension is 11.4 x 6.4cm (4.5" x 0.25 inches)

Dimensions including rod 20.3 x 3.2 x 2.5 cm

Weight 3 oz.

### Small, easy to use

Mixer is ideal for mixing samples in containers too narrow for many other mixers.

### Easy to replace batteries

2 AA batteries are easily replaceable in this mixer design.

---

## Specifications

Dimensions (H x W x D):	203 mm x 32 mm x 25 mm
Electrical Specifications:	2 AA Batteries
Length:	203 mm
Power:	Hand-Held
Quantity:	1 each
Special Features:	Stirrer, Mixer
Type:	Stirrer
Weight:	0.085 kg
What's included?:	Includes: Mixer with rod, carrying case and 2 AA batteries.

---

## What's included?

Includes: Mixer with rod, carrying case and 2 AA batteries.