



## Centrifuge, Laboratory, 230 Vac, 50 Hz

Product #:

2676502

 **Hazardous**

### **OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Versatile benchtop unit for sample preparation such as settleometer test. Stainless-steel motor shaft and permanently lubricated ball bearings provide long, trouble-free service life. Centrifuge accepts tubes with dimensions of 13 - 17 mm diameter x 100 - 133 mm long.

### **Versatile**

Can accommodate 7 mL, 10 mL, or 15 mL vacutainer tubes.

### **Built-in timer**

Allows runs up to 15 minutes.

### **Fast start**

Full speed is achieved in less than 15 seconds with a full load.

### **Multiple samples**

Six-place, fixed angle rotor with buckets and numbered positions standard in centrifuge.

### **Trouble-free service life**

Stainless steel motor shaft and permanently lubricated ball bearings provide a long operational life.

---

## Specifications

Electrical Specifications:	230 Vac, 50 Hz, 1.05 A
Power:	230 Vac, 50 Hz
Quantity:	1 each
Special Features:	Centrifuge, Laboratory
What's included?:	Includes: Centrifuge, rotor, shield and instruction manual.

---

## What's included?

Includes: Centrifuge, rotor, shield and instruction manual.