



## Pump, Vacuum, Hand-operated

<b>Product #:</b>	<b>1428300</b>
ZAR Price:	Contact Hach
Ships within 1 week	

This unit is ideal for field use and is a real convenience in the lab. Pumping rate is 15 cc (0.50 oz) per stroke. Entire unit is corrosion-resistant, sealed, and self-lubricating. Nozzle fits standard 6 mm (1/4") I.D. tubing. Supplied with 61 cm (2') of clear plastic tubing. Choose with or without vacuum gauge.

---

### Specifications

Quantity:	1 each
-----------	--------