



Branson Ultrasonic Bath, 115 Vac, 60 Hz

Product #:	2489500
ZAR Price:	Contact Hach
Call for ship date	

Cleans small and medium sized parts, including spectrophotometer cells. Unit

comes with heater, timer, and temperature control plus a 2.8 L (0.75-gal) capacity tank. Tank is resistant to most cleaning solutions.

Dimensions: exterior (L x W x D) 33.8 x 30 x 29.5 cm (13.3 x 12 x 11.6

Specifications

Electrical Specifications:	115 Vac, 60 Hz
Length:	338 mm
Quantity:	1 each