



## MultiTest kit case (145 x 206 x 66 mm), blue polypropylene

<b>Product #:</b>	<b>4660100</b>
ZAR Price:	Contact Hach
Available	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

**Recessed handle in lid provides easy stacking**

**Contents are carried in upright position**

**Includes recessed flat-panel areas for labeling**

---

### Specifications

Case Style:	B
Dimensions (H x W x D):	145 mm x 206 mm x 66 mm
Height:	5.75 in