



Iron FerroVer Color Disc, Low Range, Indoor Light, 0 - 1 mg/L

Product #:

9263600

ZAR Price:

Contact Hach

Ships within 3 days

Color Disc for measuring Iron FerroVer, 0 - 1 mg/L, with indoor light source. Used in test kit Model IR-18A (Prod. No. 146500).

The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons

Specifications

Method:	Color disc - 1, 10 Phenanthroline
Parameter:	Iron, medium range - as Fe
Range 2:	0.02 - 1.0 mg/L Fe