



Surface Water Test Kit

Product #:

2559833

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

By providing economical, easy-to-use tests for ammonia (0 - 2.4 mg/L), chlorine (0 - 3.4 mg/L), pH (0 - 14), nitrate (0 - 40 mg/L), dissolved oxygen (0.2 - 20 mg/L), phosphorus (0 - 40 mg/L), and temperature, this kit is an ideal tool for source water monitoring and wastewater treatment plants. The rugged, lightweight kit includes the convenient pH Pocket Pro Tester and all necessary reagents and apparatus for testing.

Simple, accurate methods

Utilizes drop count titration, color disc colorimetry, and pH Pocket Pal Tester

Specifications

Application: Surface Water

Platform: Color Disc

What's included?: Includes ten reagents, deionized water, pH tester, measuring tube, bottles, comparator box, BOD bottle, goggles, tubes, color discs, instruction sheet, and carrying case.

What's included?

Includes ten reagents, deionized water, pH tester, measuring tube, bottles, comparator box, BOD bottle, goggles, tubes, color discs, instruction sheet, and carrying case.