



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Product #:

223002

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

Specifications

Application:	Water Quality
Model:	CN-39WR
Number of tests:	100
Parameter:	Chlorine, Hardness, Iron, and pH
Platform:	Color Disc
Range:	0 - 3.4 mg/L