



## Test kit, boiler feed & scale, model BSC-1, drop count titration

Product #:

2350300

ZAR Price:

Contact Hach

Ships within 2 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. An economical test kit for monitoring scale formation in small industrial boilers, the BSC-1 Boiler Feed and Scale Test Kit uses drop count titration methods to determine levels of alkalinity, chloride, and hardness.

---

### Specifications

Application:	Boiler Water
Model:	BSC-1
Platform:	Drop Count Titrator