



Test kit, professional boiler treatment/boiler feed, model BT-DT

Product #:

2350600

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. This kit provides testing capability for alkalinity, chloride, hardness, sulfite, phosphate, and pH, with a wide choice of ranges.

Offers testing of treatment chemicals and feedwater parameters, with the accuracy demanded by water treatment professionals

The included Hach Digital Titrator simplifies titration procedures, saves time, and provides greater accuracy than drop count testing

Kit also includes color disc for simple, fast analysis

Specifications

Application:	Boiler Water
Platform:	Digital Titrator