



Nitrification Control Test Kit

Product #: 2882200
ZAR Price: Contact Hach
Call for ship date



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The Nitrification Control Kit is used for the routine monitoring of water quality parameters that contribute to or indicate nitrification is occurring in chloraminated water systems. This portable and versatile kit includes a DR 900 Colorimeter and reagents for Free Ammonia and Monochloramine, Nitrite, and Total Chlorine tests.

Used for the routine monitoring of water quality parameters that contribute to or indicate nitrification is occurring in chloraminated water systems

Versatile kit contains the DR 900 Colorimeter and reagents to perform tests for Free Ammonia, Monochloramine, Nitrite, Total Chlorine (low range and high range)

Specifications

Application:	Nitrification
Instrument:	Kit includes DR 900 + Sample Cells
Range:	0.02 - 0.50 mg/L $\text{NH}_3\text{-N}$
Range 2:	0.04 - 4.50 mg/L Cl_2
Range 3:	0 - 0.350 mg/L $\text{NO}_2\text{-N}$
Range 4:	0 - 2,00 mg/L Cl_2
Range 5:	0.1 - 10.0 mg/L Cl_2
What's included?:	Includes DR 900 colorimeter, five reagents, scissors, CD, sample cell and cap.

What's included?

Includes DR 900 colorimeter, five reagents, scissors, CD, sample cell and cap.