



Cyanuric acid test kit, model CY-3

Product #: 185102
ZAR Price: Contact Hach
Ships within 2 weeks

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Application:	Pool and Spa
Case Style:	D
Method:	Turbidimetry
Model:	CY-3
Number of tests:	100
Parameter:	Cyanuric Acid
Range:	20 - 100 mg/L
Smallest increment:	20
Test Kit Model:	CY-3