



## Hardness, iron, and pH test kit, model HA-62

**Product #:**

**183700**

**ZAR Price:**

**Contact Hach**

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Model HA-62

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

**Popular tests at a great price in Hach's convenient and protective test kit carrying case.**

**Uses color disc and drop count titration methods**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Application:	Water Quality
Model:	HA-62
Number of tests:	100
Parameter:	Hardness, Iron, and pH
Platform:	Color Disc
Range:	6.6 - 8.4 pH units