



Hardness, iron, and pH test kit, model HA-62A

Product #:

183701

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Application:	Water Quality
Model:	HA-62A
Number of tests:	100
Parameter:	Hardness, Iron, and pH (6.6-8.4)
Platform:	Color Disc
Range:	6.6 - 8.4 pH units