



Hardness, iron, and pH test kit, model HA-62B

Product #:

183702

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Model HA-62B

Includes color disc and color comparator for colorimetric comparison of iron and pH, as well as a drop-count titrator for hardness determination. Test kit produces a wide range of pH tests.

Combines frequently measured parameters in one rugged kit

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Application:	Water Quality
Model:	HA-62B
Number of tests:	100
Parameter:	Hardness, Iron, and pH (4-10)
Platform:	Color Disc