



Iron and manganese color disc test kit, model IR-20

Product #:

146300

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3.

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Application:	Water Quality
Model:	IR-20
Parameter:	Iron and Manganese
Platform:	Color Disc
Range:	0 - 3.0
Range 2:	0 - 4.0 mg/L Fe
Smallest Increments Steps:	0.2