



## Test kit, hardness & iron color disc, model HA-77

**Product #:**

**202300**

**ZAR Price:**

Contact Hach

**Call for ship date**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Hardness & Iron Water Quality Test Kit is the kit to use if you need better accuracy than that provided by a test cube. Utilizing drop count titration and color disc colorimetry, this kit contains 100 tests.

**The kit to use for more resolution in a visual test**

**Uses simple, economical drop count titration and color disc methods**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Application:	Water Quality
Model:	HA-77
Parameter:	Hardness and Iron
Platform:	Color Disc