



Storm Water Test Kit, Model SW-1

Product #:

2481300

ZAR Price:

Contact Hach

Call for ship date

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The Storm Water Multi-Parameter Water Quality Test Kit is ideal for on-site or laboratory monitoring of storm sewer outflow and industrial discharge. The kit includes all necessary reagents and apparatus for monitoring detergents (0-1.2 mg/L), total chlorine (0-3.4 mg/L), copper (0-4 mg/L), and phenols (0-4 mg/L), among other parameters. Also included is a pH Pocket Pro Tester for quick pH measurements.

Tailored specifically for storm water testing needs

Designed for on-site or laboratory monitoring of storm sewer outflow or industrial discharge

Useful for post-event testing

Use to screen for pollutants after rainfall or snow events to meet stormwater discharge monitoring requirements

Simple, accurate methods

Uses drop count titration and color disc methods

Specifications

Application:	Stormwater
Platform:	Color Disc
What's included?:	Includes ten reagents, demineralizer bottle, clippers, test tubes, beaker, graduated cylinder, flask, funnel, filter paper, color discs, comparator box, viewing tube, pH tester, instructions and carrying case

What's included?

Includes ten reagents, demineralizer bottle, clippers, test tubes, beaker, graduated cylinder, flask, funnel, filter paper, color discs, comparator box, viewing tube, pH tester, instructions and carrying case