



Test kit, iron & pH color disc , model IR-23

Product #:

256200

ZAR Price:

Contact Hach

Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Utilizing the simple, accurate color disc method, this Iron & pH Water Quality Test Kit provides accurate, on-site measurement of two important parameters to help maintain water quality within the desired range. Range: 0-10 mg/L/ 4-10 pH unit.

Uses color disc to provide accurate, on-site measurement of two important parameters to help maintain water quality within the desired range

Specifications

Application:	Water Quality
Model:	IR-23
Parameter:	Iron and pH
Platform:	Color Disc