



Deluxe water conditioning demonstration test kit, model AD-16E

Product #:

148410

ZAR Price:
Available

Contact Hach



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Designed to boost your testing capabilities, Hach's most complete water conditioning test kit includes tests and demonstrations for seven critical water quality parameters. A optional softener regeneration test kit is available (P/N - 2666000).

Uses color disc and drop count titration methods

In addition to soap consumption and softener demonstrations, includes tests to measure iron, pH, and nitrate

Also includes the Total Dissolved Solids (TDS) Pocket Pal™ Tester

Comes complete with a sturdy carrying case and features viewing tube, water softener, and a compartmentalized reagent cabinet to keep chemistries neat and well-organized

Specifications

Application:	Demonstration
Model:	AD-16E
Platform:	Color Disc