



## Boron High Range Starter Set

**Product #:** 2527900  
**ZAR Price:** Contact Hach  
**Ships within 1 week**



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Starter Set for determination of high range boron by the Carmine method.

Hach Method 10252. 2 - 50 mg/L Boron

Set contains BoroVer 3 Reagent Powder Pillows, Sulfuric Acid, Deionized Water, plus an Erlenmeyer flask, 25 reaction tubes and a cooling rack.

**Colorimetric method for the determination of Boron.**

---

### Specifications

Application:	Boron Testing
Range 2:	2 - 50 mg/L B
What's included?:	Includes two reagents, deionized water, flask, test tube rack, and 25 vials.

---

### What's included?

Includes two reagents, deionized water, flask, test tube rack, and 25 vials.