



## Test kit, phosphorus, orthophosphate (reactive) , model PO-19A

Product #:

224801

ZAR Price:

Contact Hach

Ships within 3 days



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Simple, low-cost color disc method is the most accurate visual test kit method

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifications

Case Style:	F
Measuring range:	Triple Range
Method:	Color disc/Ascorbic acid
Model:	PO-19A
Number of tests:	100
Parameter:	Phosphorus
Platform:	Color Disc
Range:	0 - 0.8 mg/L PO <sub>4</sub>
Range 2:	0 - 4 mg/L PO <sub>4</sub>
Range 3:	0 - 40 mg/L PO <sub>4</sub>
Smallest increment:	0.02
Standard range:	0 - 0.8 mg/L PO <sub>4</sub> ; 0 - 4 mg/L PO <sub>4</sub> ; 0 - 40 mg/L PO <sub>4</sub>
Test Kit Model:	PO-19A

(Includes filtration for turbid samples.)